

GV 350[®]

High Volt Pulsed Stimulator



"Live and Enjoy Life Pain Free"



GV 350®



PACKAGE INCLUDES

- 1 GV350®
- 4 2"x 2" electrodes
- 1 4"x 7" electrode
- 3 Lead wires
- 4 AA alkaline batteries
- 1 Instruction manual
- 1 Carrying pouch

SPECIAL FEATURES

- Timing options
- Low battery status indicator
- Polarity option switch
- Adjustable settings
- Retractable dials for safety
- Belt clip for portability
- HCPCS E0745 or E1399

GENERAL SPECIFICATIONS

SIZE	3.8" x 2.60" x 1"
WEIGHT	3.1 oz.
POWER	4 AA batteries (optional wall adapter available purchased separately)
CHANNELS	2 Active, 1 Dispersive return
MODES	Continuous and Alternating
WAVEFORM	Pulsed Twin Peak Monophasic
PULSE RATE	1 - 100 Hz (adjustable)
PULSE WIDTH	Interpulse intervals of 100 microseconds. 5 microseconds at 50% of peak voltage.
INTENSITY	0 - 700 mA peak (adjustable)
ON TIME	Continuous on <i>Continuous</i> mode or 5 seconds fixed on the <i>Alternating</i> mode.
OFF TIME	5 seconds fixed on <i>Alternating</i> mode and continuous on <i>Continuous</i> mode.
VOLTAGE	0 - 350 volts (adjustable)
POLARITY	Positive and negative active electrodes
TIMER	C, 15, 30, 60, minutes
TOLERANCES	+/- 10%

Data recorded across a 500 OHM load resistance



1954 Kellogg Avenue
 Carlsbad, CA 92008
 Ph: (800) 726-8367 USA
 Ph: (760) 579-0801 INTL
 Fx: (760) 929-9953
www.bmls.com
information@bmls.com



Cage Code 3R8H6



Optional wall adapter sold separately

BioMedical Life Systems has been producing the best portable therapy devices and accessories for over thirty-five years. The result is a comprehensive family of products meeting the needs of physical therapists, pain clinics, medical insurance, veterinarians, professional athletics, and individuals the world over.

This is the ONLY hand-held High Volt Pulsed Stimulator delivering an intensity level of 700 mA on the market today. Used by universities, professional sport organizations and hospitals worldwide. The GV 350® is a two channel device that delivers full powered twin peak alternating or continuous stimulation. It does this by delivering pulsed twin peak monophasic waves delivered in pairs for maximum performance. Incorporated with a timer, an adjustable pulse rate control, two active outputs, and one dispersive return.

Widely used for it's wound healing and pain management results as well as other indications for edema and muscle strengthening.

This unit is powered by four AA batteries though an optional wall adapter is available to purchase separately and includes everything needed to manage your therapy sessions.