

Impulse® 4000T

Analog TENS Stimulator



BioMEDICAL
LIFE • SYSTEMS

Therapy Solutions
since 1984

"Live and Enjoy Life Pain Free"



Impulse® 4000T



PACKAGE INCLUDES

- 1 Impulse® 4000T
- 4 2"x 2" electrodes
- 2 Lead wires
- 1 9V alkaline batteries
- 1 Instruction manual
- 1 Carrying pouch

SPECIAL FEATURES

- Three modes of operation
- Compact light weight design
- Safety dial cover
- Fully adjustable settings
- Belt clip
- HCPCS E0730

GENERAL SPECIFICATIONS

SIZE:	3.8" x 2.60" x 1"
WEIGHT:	3.35 oz
POWER:	9V Battery, E-block, type 6F22
CHANNELS:	Dual (2)
WAVEFORM:	Asymmetrical, biphasic square-wave
PULSE RATE:	2 - 150 Hz (Hertz or pps) adjustable
PULSE WIDTH:	50 - 250 microseconds (µS) adjustable
INTENSITY:	0- 98 mA peak to peak
OUTPUT:	Constant current
MODES:	<u>Continuous</u> - Pulse Rate/Pulse Width are adjustable. <u>Burst</u> - 8 Pulses per burst, 2 bursts per second. Pulse Width is adjustable. <u>Modulation</u> - Decrease/increase of set width of 50% over a 5-second period.
TIMER:	Continuous, 15, 30, 60 minutes
TOLERANCES:	+/- 10%

Data recorded across a 500 OHM load resistance

Does pain management have to involve drugs? Absolutely not! Pain can affect your mood, social events, and joyful activities with family & friends. Researchers agree that the use of an electrotherapy device can provide sufficient relief from neuropathy, muscle, or joint pain so you can live and enjoy life!

The Impulse® 4000T is a portable dual channel transcutaneous electrical nerve stimulator (TENS). It utilizes small electrical currents affecting the nerves from sending pain signals to the brain so that pain is not perceived. Indicated for chronic pain, post surgical pain, acute pain, and neuropathy relief.

This TENS unit is ergonomically designed and easy to use. Incorporate three modes of operation; constant, burst, and modulated. Parameters for each mode are fully adjustable. This TENS unit with timer is perfect for those seeking timed sessions.

This design incorporates a flip top amplitude control cover for safety keeping the intensity from changing accidentally. The 9 volt battery compartment has been designed for simple insertion and removal of the battery. To provide comfort while wearing, the belt clip is provided.

This device comes with all required accessories needed to manage your therapy sessions.



1954 Kellogg Avenue
Carlsbad, CA 92008
Ph: (800) 726-8367 USA
Ph: (760) 579-0801 INTL
Fx: (760) 929-9953
www.bmls.com
information@bmls.com



Cage Code 3R8H6

